



LC-MS/MS Assays

Product Overview

Complementary Medicine

Gastroenterology / Nutrition

Skeletal System / Bone Metabolism

Clinical Immunology

TDM





Product *	Cat. No.	Matrix	Volume	Mode	Polarity
1,25-(OH) ₂ -Vitamin D ₃ /D ₂ ImmuTube® LC-MS/MS Kit	KM1000	Plasma, Serum	500 µl	MRM	ESI positive
1,25-(OH) ₂ -Vitamin D ₃ /D ₂ IDKenhance® ImmuTube® LC-MS/MS Kit	KM1002	Plasma, Serum	500 µl	MRM	ESI positive
1,25-(OH) ₂ -Vitamin D ₃ /D ₂ ImmuPlate® DIY LC-MS/MS Kit	KM1020	Plasma, Serum	500 µl	MRM	ESI positive
1,25-(OH) ₂ -Vitamin D ₃ /D ₂ IDKenhance® ImmuPlate® DIY LC-MS/MS Kit	KM1022	Plasma, Serum	500 µl	MRM	ESI positive
1,25-(OH) ₂ -Vitamin D ₃ /D ₂ IDKmagnetics® DIY LC-MS/MS Kit **	KM1040	Plasma, Serum	500 µl	MRM	ESI positive
1,25-(OH) ₂ -Vitamin D ₃ /D ₂ IDKenhance® IDKmagnetics® DIY LC-MS/MS Kit **	KM1042	Plasma, Serum	500 µl	MRM	ESI positive
1,25-(OH) ₂ -Vitamin D ₃ /D ₂ ImmuTube® LC-MS/MS Extraction Kit	KM1100	Plasma, Serum			<i>for 50 samples</i>
1,25-(OH) ₂ -Vitamin D ₃ /D ₂ IDKenhance® ImmuTube® LC-MS/MS Extraction Kit	KM1102	Plasma, Serum			<i>for 50 samples</i>
1,25-(OH) ₂ -Vitamin D ₃ /D ₂ ImmuPlate® DIY LC-MS/MS Extraction Kit	KM1120	Plasma, Serum			<i>for 96 samples</i>
1,25-(OH) ₂ -Vitamin D ₃ /D ₂ IDKenhance® ImmuPlate® DIY LC-MS/MS Extraction Kit	KM1122	Plasma, Serum			<i>for 96 samples</i>
1,25-(OH) ₂ -Vitamin D ₃ /D ₂ IDKmagnetics® DIY LC-MS/MS Extraction Kit **	KM1140	Plasma, Serum			<i>for 96 samples</i>
1,25-(OH) ₂ -Vitamin D ₃ /D ₂ IDKenhance® IDKmagnetics® DIY LC-MS/MS Extraction Kit **	KM1142	Plasma, Serum			<i>for 96 samples</i>
ImmuTube® Vitamin D duo LC-MS/MS Kit for the determination of 24,25-(OH) ₂ -D ₃ and 25-OH-D ₃ /D ₂	KM1310	Plasma, Serum	50 µl	MRM	ESI positive
ImmuTube® Vitamin D duo LC-MS/MS Extraction Kit	KM1311	Plasma, Serum			<i>for 50 samples</i>
IDKmagnetics® Vitamin D duo LC-MS/MS Kit for the determination of 24,25-(OH) ₂ -D ₃ and 25-OH-D ₃ /D ₂	KM1320	Plasma, Serum	50 µl	MRM	ESI positive
IDKmagnetics® Vitamin D duo LC-MS/MS Extraction Kit	KM1321	Plasma, Serum			<i>for 96 samples</i>
Pregnenolone sulfate LC-MS/MS Kit	KM2000	Plasma, Serum	50 µl	MRM	ESI negative



Product *	Cat. No.	Matrix	Volume	Mode	Polarity
IDK® Cortisol & Cortisone LC-MS/MS Kit	KM2110	Urine	300 µl	MRM	ESI positive
IDK® Metanephrines LC-MS/MS Kit	KM2200	EDTA Plasma	500 µl	MRM	ESI positive
IDK® Methylmalonic acid LC-MS/MS Kit	KM3110	Plasma, Serum, Urine	100 µl	MRM	ESI negative
3-OH-Isovalerylcarnitine LC-MS/MS Kit	KM3200	Citrate Plasma	50 µl	MRM	ESI positive
Hepcidin-25 LC-MS/MS Kit	KM4000	Plasma, Serum	200 µl	MRM	ESI positive
IDK® Vitamin B₁ & B₆ LC-MS/MS Kit	KM5100	Whole blood	50 µl	MRM	ESI positive
IDK® Vitamin A & E LC-MS/MS Kit	KM5200	Plasma, Serum	25 µl	MRM	ESI positive
IDK® Psychotropic Medication 1 LC-MS/MS Kit 11 analytes: Clozapine, Norclozapine, Nortriptyline, E-10-OH-Nortriptyline, Z-10-OH-Nortriptyline, Amitriptyline, Clomipramine, Norclomipramine, Imipramine, Desipramine, Zuclopenthixol	KM6310	Serum	25 µl	MRM	ESI positive
IDK® Casomorphin/Gliadorphin Peptides LC-MS/MS Kit	KM8000	Stabilized urine	1 ml	MRM	ESI positive
IDK® Immunosuppressants LC-MS/MS Kit Tacrolimus, Cyclosporin A, Sirolimus, Everolimus	KM9100	Whole blood	50 µl	MRM	ESI positive
IDKmonitor® Vedolizumab drug level LC-MS/MS Kit (e.g. Entyvio®) <i>For research only!</i>	KMR9600	Serum	10 µl	MRM	ESI positive

** coming soon

Continuous external quality control of our LC-MS/MS assays in interlaboratory tests:

INSTAND e.V., RfB and DEQAS

Our assays are calibrated and produced with reference materials.

(NIST, CDC, Cerilliant and more)

* For laboratory professional use only.

All products are for research use only in the U.S. Not for use in diagnostic procedures.

Contact your sales representative for more information.

Contact our LC-MS team!



Our service & support:

- Fast feedback to all customers by e-mail, phone, remote or on site
- Application notes for different LC-MS systems are provided with the kit
- Product presentation, demonstration and training on site

LCMS@immundiagnostik.com
+49 6251 70190-242